AT51004

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USER MANUAL ROLLATOR

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MANUEL D'UTILISATION ROLLATOR

HANDLEIDING ROLLATOR

MANUAL DE USO RODILLO

MANUALE D'USO RULLO

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INTENDED USE: The product is intended to support the movement of people with reduced walking ability. The walker is designed for one person.

CONTRAINDICATIONS: Physical or mental limitations (e.g. visual impairment) that prevent the safe handling of the walker.

Unpacking the shipment

Walker with accessories (two shopping bags included)

fig.1

Unfolding the frame

- A. Unfold the frame by pressing on the side parts of the seat.
- B. Fold the frame by pulling the strap on the seat.

fig.2

The height of the handles is adjusted with the help of knobs

To turn the knob and adjust the height, make sure the knob is tightened before pushing the handle.


fig.3

ASSEMBLING

All the steps described above should be carried out in reverse order.

CAUTION: When using and operating the walker, and when folding it and adjusting the mechanisms, there may be a risk of entrapment and / or pinching of body parts of the user / accompanying person in opening / in gaps between the components. These steps should be performed with extreme caution. After completing the adjustment, stabilize the position by tightening the nuts / bolts carefully.

ATTENTION:

In the event of pain, allergic reactions or other disturbing symptoms, unclear for the user, related to the use of a medical device, consult a healthcare professional.

ATTENTION:

When a product-related „serious incident“ occurs, which directly or indirectly led, or could lead to, any of the following events:

- a) death of a patient, user or other person or
 - b) temporary or permanent deterioration in the health of a patient, user or other person, or
 - c) a serious threat to public health
- the above „serious incident“ must be reported to the manufacturer and to the competent authority of the Member State in which the user or patient is resident. In the case of Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

ATTENTION: The manufacturer is not responsible for damage caused by neglect of maintenance, inadequate servicing, or resulting from failure to comply with the recommendations contained in this manual.

ATTENTION: it is forbidden to use the product other than in accordance with its intended purpose!

CAUTION: There may be a risk of tipping over if the product is used improperly. Please follow the sitting / standing / moving instructions. After regulating the product stabilize the position by tightening the nuts / bolts carefully.

Braking

A. Reduce speed or stop the rollator

Pull the handles up until they stop or slow down the rollator

B. Blocking the wheels

Stand between the handles and press down on the rounded part with your hand. The wheels will lock when you hear a loud noise clicks.

C. When the walking frame is stationary, the brakes must be locked. Stand between the handles and press the bottom part of the loop with your hand, moving it back to its original position.

D. Releasing the wheel lock.

fig. 4

fig. 5

fig. 6

Technical specifications

Height: 72 - 97 cm

8" diameter wheels

Seat height: 53 cm

Distance between handles: 33 cm

Distance between rear wheels: 58 cm

Distance between rear wheel and front wheel:
49 cm

Weight of walker: 6 kg

Maximum bag load: 5 kg

Maximum load of walker: 136 kg

THIS MARK INDICATES THE
MAXIMUM USER WEIGHT.

CLEANING

Clean the product with a damp cloth and mild detergent.

MAINTENANCE

1. Please ensure that all attachments are always secured.
2. Check all mounting bolts, tighten if necessary.
3. Replace any damaged or destroyed parts immediately.

WHEEL MAINTENANCE

1. Wheels should be cleaned with warm water and a mild detergent. Pat dry with a dry cloth.
2. Lubricate periodically with 3-in-1 grease if the wheels begin to squeak.

TARGET PATIENT GROUP

Persons suffering from diseases, dysfunctions or injuries for the treatment/rehabilitation or compensation of which this device is intended (see section on intended use of the device in this instruction manual). The device can be purchased by the user on their own or on the recommendation of a doctor, therapist or other specialist. Whether you purchase the device yourself or on the recommendation of a doctor/therapist/other specialist, you must take into account the available sizes/necessary functions and variants of the device, the indications and contraindications for use and the information provided by the manufacturer.

DISPOSAL OF THE PRODUCT

The product must be disposed of in accordance with the applicable disposal conditions.

STORAGE

Store the product in a dry cool room, away from moisture and direct sunlight.

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WARRANTY CARD

Model:
Serial number:

Point of sale stamp
and a legible signature of the Seller:

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawisłańska 43, hereby grants a 12-month warranty on the rollator, starting from the date of purchase by the Buyer. Products financed by the National Health Fund are covered with the extended warranty period, which covers half of the useful life specified in the Regulation on Guaranteed Services, provided that the useful life is specified in the Regulation in months or years.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, send it immediately to ANTAR (as specified in point 1) or contact the point of sale.
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum permissible weight of the user as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts subject to normal wear and tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED. IN ORDER TO EXTEND WARRANTY PERIOD FOR PRODUCTS FINANCED BY THE NATIONAL HEALTH FUND, ATTACH A COPY OF THE ORDER FOR THE SUPPLY OF MEDICAL DEVICES BEING ORTHOPAEDIC ITEMS AND AIDS.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR Sp. J.
03-068 Warsaw, ul. Zawisłańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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